

(3) 取水口										
採取日	核種濃度 (Bq/L ※PuはmBq/L)									備考
	全β放射能	¹³⁴ Cs	¹³⁷ Cs	³ H			⁹⁰ Sr	²³⁸ Pu	²³⁹⁺²⁴⁰ Pu	
				減圧蒸留法	迅速分析	電解濃縮法				
H1. 6. 14	ND	ND	ND	0.66						
H1. 10. 24	0.03	ND	ND	1.5						
H2. 7. 3	ND	ND	ND	ND						
H2. 11. 2	ND	ND	ND	1.0						
H3. 6. 20	ND	ND	ND	ND						
H3. 11. 7	ND	ND	ND	ND						
H4. 6. 26	0.01	ND	ND	0.74						
H4. 10. 29	ND	ND	ND	1.0						
H5. 6. 25	0.01	ND	ND	ND						
H5. 10. 22	ND	ND	ND	ND						
H6. 6. 7	ND	ND	0.003	ND						
H6. 11. 5	0.01	ND	0.002	ND						
H7. 6. 6	0.01	ND	0.003	ND						
H7. 10. 4	0.01	ND	0.003	ND						
H8. 6. 6	0.04	ND	0.003	ND						
H8. 10. 4	0.02	ND	0.003	ND						
H9. 6. 5	ND	ND	0.002	0.60						
H9. 10. 9	0.02	ND	ND	ND						
H10. 6. 5	0.01	ND	0.002	ND						
H10. 10. 13	ND	ND	0.001	0.50						
H11. 5. 26	0.01	ND	0.002	ND						
H11. 10. 12	0.01	ND	ND	ND						
H12. 5. 30	ND	ND	0.003	1.3						
H12. 10. 10	ND	ND	0.003	ND						
H13. 5. 18	ND	ND	0.001	ND						
H13. 10. 16	0.02	ND	0.002	0.58						
H14. 5. 20	ND	ND	0.002	ND						
H14. 10. 25	0.02	ND	0.001	ND						
H15. 5. 21	ND	ND	ND	ND						
H15. 10. 20	ND	ND	0.002	ND						
H16. 5. 19	0.02	ND	ND	ND						
H16. 10. 28	0.02	ND	ND	ND						
H17. 5. 23	ND	ND	ND	ND						
H17. 10. 26	0.01	ND	ND	ND						
H18. 5. 17	0.01	ND	0.001	ND						
H18. 10. 12	0.02	ND	0.002	ND						
H19. 5. 16	0.02	ND	0.001	ND						
H19. 10. 10	0.01	ND	0.002	0.46						
H20. 5. 16	0.01	ND	0.001	0.37						
H20. 10. 17	ND	ND	ND	ND						
H21. 5. 15	0.01	ND	0.001	ND						
H21. 10. 13	ND	ND	0.001	ND						
H22. 5. 21	ND	ND	ND	ND						
H22. 10. 7	ND	ND	0.001	ND						
H23. 4. 1	欠測	欠測	欠測	欠測						
H23. 10. 1	欠測	欠測	欠測	欠測						
H25. 7. 31	0.02	ND	ND	ND			0.015	ND	ND	
H25. 8. 19	1.7	0.29	0.75	3.0			2.9	ND	ND	
H25. 9. 19	0.74	0.96	2.0	4.9			0.95	ND	ND	
H25. 10. 5	0.02	ND	0.073	0.53			0.005	ND	ND	
H25. 10. 30	1.2	0.47	1.1	6.2			1.6	ND	ND	
H25. 11. 12	0.07	0.083	0.21	0.40			0.1	ND	0.006	
H25. 12. 9	0.15	0.15	0.32	0.76			0.22	ND	ND	
H26. 1. 14	0.18	0.12	0.33	0.80			0.21	ND	ND	
H26. 2. 3	0.68	0.54	1.3	3.4			1.1	ND	0.006	
H26. 3. 10	0.07	ND	0.15	ND			0.046	ND	0.01	
H26. 4. 14	0.03	ND	ND	ND			0.003	ND	0.008	

(3) 取水口

採取日	核種濃度 (Bq/L ※PuはmBq/L)									備考
	全β放射能	¹³⁴ Cs	¹³⁷ Cs	³ H			⁹⁰ Sr	²³⁸ Pu	²³⁹⁺²⁴⁰ Pu	
				減圧蒸留法	迅速分析	電解濃縮法				
H26.5.19	0.05	ND	ND	ND			0.005	ND	ND	
H26.6.16	0.03	ND	ND	ND			0.003	ND	ND	
H26.7.7	0.38	0.24	0.72	2.6			0.66	ND	ND	
H26.8.20	0.19	0.12	0.25	1.1			0.26	ND	ND	
H26.9.16	0.23	0.35	0.94	2.5			0.4	ND	ND	
H26.10.20	0.08	0.11	0.33	0.61			0.12	ND	ND	
H26.11.10	0.17	0.13	0.36	1.7			0.14	ND	ND	
H26.12.8	0.05	ND	ND	ND			0.008	ND	ND	
H27.1.14	0.12	ND	0.082	ND			0.052	ND	ND	
H27.2.10	0.18	0.097	0.16	0.59			0.47	ND	ND	
H27.3.3	0.14	0.079	0.23	1.4			0.23	ND	ND	
H27.4.22	0.02	ND	ND	ND			0.002	ND	ND	
H27.5.18	0.02	ND	ND	ND			0.002	ND	ND	
H27.6.16	0.03	ND	ND	ND			0.008	ND	ND	
H27.7.21	0.37	0.22	0.83	2.6			0.68	ND	ND	
H27.8.17	0.19	0.13	0.48	1.1			0.46	ND	ND	
H27.9.14	0.03	ND	0.087	ND			0.002	ND	ND	
H27.10.13	0.19	0.1	0.46	2.1			0.27	ND	ND	
H27.11.17	0.02	ND	0.07	ND			0.01	ND	0.007	
H27.12.14	0.06	ND	0.081	ND			0.003	ND	0.012	
H28.1.22	0.04	ND	ND	ND			0.002	ND	0.007	
H28.2.8	0.04	ND	0.15	0.63			0.014	ND	ND	
H28.3.3	0.04	ND	0.12	ND			0.005	ND	0.01	
H28.4.20	0.03	0.003	0.013	ND			0.0011	ND	ND	
H28.5.16	0.02	0.006	0.031	ND			0.0013	ND	0.01	
H28.6.15	0.02	0.002	0.015	ND			0.0009	ND	ND	
H28.7.11	ND	ND	0.011	ND			ND	ND	ND	
H28.8.3	0.02	ND	0.01	ND			0.0005	ND	ND	
H28.9.15	0.12	0.2	1.1	0.55			0.087	ND	ND	
H28.10.18	0.03	0.005	0.032	ND			0.0006	ND	ND	
H28.11.15	0.03	0.038	0.22	0.65			0.013	ND	ND	
H28.12.12	0.03	0.021	0.13	ND			0.0056	ND	ND	
H29.1.20	0.03	0.006	0.034	ND			0.0026	ND	0.008	
H29.2.14	0.03	0.038	0.24	ND			0.0091	ND	ND	
H29.3.7	0.04	0.004	0.025	ND			0.0013	ND	0.007	
H29.4.20	0.03	0.021	0.15	ND			0.0071	ND	ND	
H29.5.16	0.03	0.011	0.08	ND			0.0067	ND	0.009	
H29.6.13	0.02	ND	0.006	ND			0.0007	ND	ND	
H29.7.10	0.05	0.016	0.12	ND			0.005	ND	ND	
H29.8.18	0.03	0.009	0.067	0.58			0.0085	ND	ND	
H29.9.14	0.02	ND	0.008	ND			0.0011	ND	ND	
H29.10.17	0.02	ND	0.012	ND			0.0011	ND	ND	
H29.11.14	0.02	0.003	0.022	ND			0.002	ND	ND	
H29.12.5	0.03	0.018	0.14	0.43			0.01	ND	ND	
H30.1.16	0.02	0.002	0.024	0.45			0.0019	ND	ND	
H30.2.13	0.03	ND	0.008	ND			0.0011	ND	0.005	
H30.3.13	0.02	0.013	0.12	ND			0.017	ND	ND	
H30.4.20	0.02	ND	0.01	ND			0.0007	ND	ND	
H30.5.16	0.02	0.008	0.086	ND			0.013	ND	ND	
H30.6.14	0.03	0.008	0.071	ND			0.01	ND	0.007	
H30.7.10	0.02	ND	0.007	ND			0.0008	ND	ND	
H30.8.19	0.03	0.003	0.045	ND			0.0012	ND	ND	
H30.9.13	0.03	0.031	0.34	0.66			0.013	ND	0.008	
H30.10.5	0.02	0.012	0.14	0.44			0.01	ND	0.003	
H30.11.14	0.02	ND	0.016	ND			0.0008	ND	0.009	
H30.12.11	0.02	0.004	0.032	ND			0.0011	ND	ND	
H31.1.17	0.03	ND	0.02	ND			0.0008	ND	ND	

(3) 取水口

採取日	核種濃度 (Bq/L ※PuはmBq/L)									備考
	全β放射能	¹³⁴ Cs	¹³⁷ Cs	³ H			⁹⁰ Sr	²³⁸ Pu	²³⁹⁺²⁴⁰ Pu	
				減圧蒸留法	迅速分析	電解濃縮法				
H31. 2. 13	0.02	ND	0.031	ND			0.0012	ND	0.007	
H31. 3. 18	0.03	ND	0.02	ND			0.0011	ND	ND	
H31. 4. 17	0.03	ND	0.032	ND			0.0012	ND	0.009	
R1. 5. 10	0.02	ND	0.006	ND			0.0006	ND	ND	
R1. 6. 4	0.03	0.006	0.066	ND			0.0026	ND	0.009	
R1. 7. 2	0.02	0.028	0.38	0.51			0.013	ND	ND	
R1. 8. 1	0.02	0.014	0.18	0.51			0.0047	ND	ND	
R1. 9. 20	0.02	0.023	0.33	0.66			0.010	ND	ND	
R1. 10. 2	0.03	0.008	0.15	0.83			0.0034	ND	ND	
R1. 11. 21	0.02	0.003	0.048	ND			0.0014	ND	ND	
R1. 12. 11	0.02	0.009	0.14	0.82			0.0045	ND	ND	
R2. 1. 8	0.05	0.012	0.17	0.89			0.0051	ND	ND	
R2. 2. 4	0.03	0.003	0.046	ND			0.0010	ND	0.010	
R2. 3. 12	0.03	ND	0.026	ND			0.0013	ND	ND	
R2. 4. 22	0.04	ND	0.021	ND			0.0009	ND	0.007	
R2. 5. 14	0.03	0.003	0.054	0.39			0.0029	ND	ND	
R2. 6. 2	0.02	0.004	0.065	ND			0.0025	ND	0.009	
R2. 7. 3	0.02	ND	0.006	ND			0.0006	ND	ND	
R2. 8. 6	0.03	ND	0.009	ND			0.0006	ND	ND	
R2. 9. 11	0.03	ND	0.005	ND			0.0010	ND	ND	
R2. 10. 20	0.02	ND	0.011	ND			0.0009	ND	ND	
R2. 11. 12	0.03	ND	0.017	ND			0.0008	ND	ND	
R2. 12. 4	0.03	ND	0.032	ND			0.0009	ND	ND	
R3. 1. 7	0.04	0.004	0.056	ND			0.0020	ND	ND	
R3. 2. 12	0.04	ND	0.030	ND			0.0011	ND	0.014	
R3. 3. 4	0.03	0.003	0.098	ND			0.0072	ND	0.018	
R3. 4. 20	0.02	0.008	0.019	0.46			0.0096	ND	0.011	
R3. 5. 12	0.02	ND	0.022	ND			0.0015	ND	ND	
R3. 6. 3	0.02	ND	0.048	ND			0.0030	ND	0.007	
R3. 7. 6	0.02	ND	0.010	ND			0.0012	ND	ND	
R3. 8. 4	0.01	0.004	0.12	0.34			0.0046	ND	0.008	
R3. 9. 2	0.03	0.010	0.31	1.4			0.035	ND	ND	
R3. 10. 15	0.02	0.006	0.12	ND			0.0076	ND	0.009	
R3. 11. 4	0.01	ND	0.067	0.46			0.0051	ND	0.009	
R3. 12. 14	0.02	ND	0.020	ND			0.0008	ND	ND	
R4. 1. 13	0.02	0.002	0.054	ND			0.0016	ND	0.007	
R4. 2. 3	0.02	0.002	0.067	ND			0.0016	ND	ND	
R4. 3. 3	0.02	ND	0.063	ND			0.0014	ND	0.009	
R4. 4. 13	0.02	ND	0.016	ND			0.0016	ND	ND	
R4. 5. 19	0.02	ND	0.16	0.62		0.66	0.012	ND	0.010	
R4. 6. 19	0.02	ND	0.013	ND			0.0019	ND	ND	
R4. 7. 5	0.02	ND	0.007	ND			0.0008	ND	ND	
R4. 8. 2	0.01	ND	0.004	ND		0.09	0.0007	ND	ND	
R4. 9. 13	0.02	ND	0.073	0.61			0.0025	ND	ND	
R4. 10. 21	0.02	ND	0.10	0.55			0.015	ND	0.009	
R4. 11. 8	0.02	ND	0.13	0.58		0.56	0.0051	ND	0.014	
R4. 12. 9	0.05	ND	0.012	ND			0.0014	ND	ND	
R5. 1. 13	0.06	ND	0.042	ND			0.0012	ND	ND	
R5. 2. 7	0.06	ND	0.010	ND		0.05	0.0011	ND	0.012	
R5. 3. 7	0.02	ND	0.070	0.51			0.0025	ND	0.007	
R5. 4. 25	0.02	ND	0.006	ND			0.0008	ND	ND	
R5. 5. 10	0.02	0.003	0.14	ND		0.21	0.012	ND	0.007	
R5. 6. 7	0.02	0.003	0.12	0.39			0.0072	ND	0.011	
R5. 7. 11	0.02	ND	0.16	ND			0.0011	ND	ND	
R5. 8. 8	0.02	ND	0.055	0.46		0.20	0.0018	ND	ND	
R5. 8. 25					ND					
R5. 8. 30					ND					

(3) 取水口										
採取日	核種濃度 (Bq/L ※PuはmBq/L)									備考
	全β放射能	¹³⁴ Cs	¹³⁷ Cs	³ H			⁹⁰ Sr	²³⁸ Pu	²³⁹⁺²⁴⁰ Pu	
				減圧蒸留法	迅速分析	電解濃縮法				
R5. 9. 3	0.01	ND	0.004		ND	0.09	ND	ND	ND	
R5. 9. 12					ND					
R5. 9. 19					ND					
R5. 9. 26					ND					
R5. 10. 8					ND					
R5. 10. 12	0.02	0.003	0.12		ND	0.25	0.0032	ND	ND	
R5. 10. 20					ND					
R5. 10. 24					ND					
R5. 11. 3					ND					
R5. 11. 9	0.02	ND	0.072		ND	0.47	0.0030	ND	ND	
R5. 11. 14					ND					
R5. 11. 22					ND					
R5. 11. 28					ND					
R5. 12. 5	0.02	ND	0.019		ND	0.14	0.0010	ND	ND	
R5. 12. 15					ND					
R5. 12. 20					ND					
R6. 1. 10					ND					
R6. 1. 18	0.02	ND	0.012		ND	0.09	0.0015	ND	ND	
R6. 1. 31					ND					
R6. 2. 9	分析中	分析中	分析中		ND	分析中	分析中	分析中	分析中	
R6. 2. 15					ND					
R6. 3. 15	分析中	分析中	分析中		ND	分析中	分析中	分析中	分析中	